

Willaid Health Care Equipment

Sales, Service, Hire



Mobility Scooter User Checklist

- Ensure seating, tiller adjustment, mirrors and back rest are in a comfortable position.
- ALWAYS board scooter with key in the off position.
- In the event of a breakdown, disengage the drive and move the scooter onto a trailer or other vehicle
 - *Do not attempt to push for long distance.*
- Charge regularly, you CAN NOT overcharge
- Store in a dry area, wipe down if caught in the rain.
- Keep tyre pressure to recommended level
- Annual servicing can prevent scooter problems
- Read user guide for further help



Willaid Health Care Equipment Service
Hotline

02 6651 7900